



Custom Wire technologies, Inc.

INNOVATORS IN MICROFINE WIRE FOR MEDICAL INSTRUMENTATION

ISO 13485 : 2016

Custom Wire technologies offers exceptional R&D services

Custom Coils in a Week

QwikCoil



UNIFILAR COILS ONLY

- Coil ID from 0.006" and up
- Lengths up to 60"
- Closed wound or open pitched available
- Standard end cuts (axial)

PRICING

- 304V SS - \$1,250 (20 coils)
- Nitinol #1 - \$2,000 (20 coils)

*Material must be in-stock to qualify.
Customer supplied material is acceptable.
Certification included.*

STOCKED MATERIALS

304V SS (spring temper) & Nitinol #1 (straight annealed)

| 304V SS & Nitinol Round Wire | 304V SS Flat Wire | Nitinol Flat Wire |
|------------------------------|---|---|
| .001" to .010" | .001" x .005" .002" x .005" .002" x .008" .002" x .012" .003" x .009" .004" x .012" .006" x .015" | .001" x .002" .001" x .003" .002" x .005" .002" x .008" .003" x .010" |

Inquire about additional size or material availability.

For all other orders please call CWT or visit our website.

CABILITIES Overview

Medical Wire Coiling & Winding

Precision Wire Grinding

Nitinol Coils & Shape Setting

Catheter Reinforcement Coils

Medical Compression Springs

Laser & Plasma Welding

Cleanroom Assembly Services

Laser Marking & Ink Printing Services



**Custom Wire
technologies, Inc.**

INNOVATORS IN MICROFINE WIRE FOR MEDICAL INSTRUMENTATION

ISO 13485 : 2016

262-268-9388 • CUSTOMWIRETECH.COM

1123 MINERAL SPRINGS DRIVE • PORT WASHINGTON, WI 53074